

DISPO GEL

Reusable & Completely Atoxic - Hot & Cold Gel Pack



YOUR DISPOSABLE EXCELLENCE

www.dispotech.com



1 min



5 min



2 hr



REUSABLE
HOT/COLD GEL PACK
DISPO GEL

MADE IN
Italy



CLASS I MEDICAL DEVICE

* For external use only

❖ HOT INDICATION :

Lumbago, Stiff neck, Rheumatism, Hypothermia and whenever hot application is required.

❖ COLD INDICATION :

Haematomas, Sprains, Injuries, Hyperthermia and whenever cold application is required.

INSTRUCTIONS FOR USE - HOT :

HOT WATER : Immerse the gelatine bag in the boiling water for max 5 mins. Take out, wipe it and then place it into the pouch or wrap it up in a cloth and apply when required on the affected area.

MICROWAVE OVEN : Put the gelatine bag in the microwave oven for max 1 min at 650 W. Take it out and place it into the pouch or wrap it up in a cloth and apply when required on the affected area.

INSTRUCTIONS FOR USE - COLD :

FREEZER : Put the gelatine bag in the freezer for at least 2 hrs. Take it out & put it in the pouch or wrap it up in a cloth and apply when required on the affected area.

	Without cover	With cover & show box
14 X 24 CM	Rs. 488 /- pc	Rs. 1070 /- pc
14 X 18 CM	Rs. 378 /- pc	Rs. 977 /- pc
	25 pcs	1 pc

EASY ICE

Instant – Disposable ICE Pack

EASY ICE is an instant & disposable ice pack, ideal for medical, sports and domestic application which relieves pain and swelling in few minutes.

On pressing the water bubble in the middle of the pack, it reaches a temperature of -4°C (under laboratory test at 20°C)

INDICATION : Suitable in case of **Haematomas, Sprains, Injuries, Hyperthermia** and whenever cold application is required.

INSTRUCTIONS FOR USE : Break the **inside water bubble** by squeezing the bag with both hands. Then shake it in order to mix the contents. Apply on the affected or injured area.

14 X 18 CM	Without show box Rs. 131 /- pc
	25 pcs



CLASS IIA MEDICAL DEVICE

* For external use only
One time use only



MADE IN
Italy



ONE TIME USE - COLD PACK
EASY ICE

CAUTION / PRECAUTION

- For **Dispogel** - Keep the gelatine bag always in its pouch or wrap it in a cloth to avoid burn • For **Easy ICE** - Temperature reaches -4°C (under laboratory test at 20°C) • For further advice please consult your doctor / physician
- Do not swallow the contents and tear the gel bag • Keep out of reach of children

Sole Importer & Distributor in INDIA



PUSHPANJALI medi India Pvt Ltd

MEDICAL COMPRESSION | ORTHOPEDICS | PHYSIOTHERAPY | REHABILITATION



We serve all!
An ISO 9001:2015 Company



SCAN FOR MORE DETAILS