



# Butler & Accessories

Donning/Doffing aid for stockings and armsleeves

medi



Steel Rs 1600/-pc  
Indian made white color Iron Rs 1275/-pc

**Long-Handle Butler**  
(for self stocking application)



Indian made Butler Short Hand-Rs.950/-pc

**Short-Handle Butler**  
(for Armsleeves)



Indian made white color Iron Rs 2100/-pc

**Hospital Butler**  
(for Clinical use)



medi rubber Gloves

Rubber Gloves Rs 575/-pc



medi Easy On Slipper

Rs 575/-pc

## SIMPLY PUT ON COMPRESSION STOCKINGS



1

Lay the compression stocking in the stirrup with the toe of the stocking facing forwards. Turn the upper border of the stocking inside out and over the frame



2

Slide the stocking down over the stirrup until the heel appears at the edge of the stirrup



3

Now place the medi Butler on the floor. Slip your foot, toes first, into the toe of the stocking. It is recommended that you lean against something firm (table, wall) for stabilisation



4

Now step into the rigid stand. Please check that the toe and the heel are correctly positioned



5

By alternately pulling and releasing the two handles you can now slide the stocking up towards your knee



6

When the stocking is just below the knee tilt the medi butler slightly backwards and remove it from the stocking



7

To finish, smooth the compression stocking upwards evenly. Please make sure you ease the stocking up bit by bit without gathering it up or pulling on the top



**PUSHPANJALI** medi India Pvt Ltd

MEDICAL COMPRESSION | ORTHOPEDICS | PHYSIOTHERAPY | REHABILITATION



SCAN FOR MORE DETAILS