M.4s® PCL dynamic



Rigid orthosis for functional posterior cruciate ligament therapy







BACK VIEW



PRODUCT BENEFITS

1 PCL Protection Support Pad

The pad presses actively on the tibia (shinbone) and holds it in the desired position.



2 Individually adjustable tensioning dial

The integrated tensioning dial enables the support pad to be adjusted exactly to suit the individual patient.



The flat construction of the hinge enables the orthosis to be worn under clothing.



4) Tool-free wedge change

Physioglide TF: the joint construction with tweezers in each hinge cover provides for tool-free wedge changes

5 Pre-formed tibia pad

The pad distributes the pressure comfortably.



I	f o	ffers : High co	mfort in wear in turn contributes to	patient compliance - Thus to Successful therapy.

.5	circumference	thigh	31 – 37	37 – 45	45 – 51	51 – 57	57 – 65	65 – 74
m 6	(cm) size		XS	S	М	L	XL	XXL

Limitations: Extension - 0°,10°,20°,30°,45° |

Flexion: 0°,10°,20°,30°,45°,60°,75°,90° Immobilisation: 0°,10°,20°,30°,45°



* PUSHPANJALI medi India Pvt Ltd 2

We serve all!



QUALITY CERTIFICATES l prices are subject to change without notice, and are not guranteed. Please check a