

The economical solution for venous problems

QUALITY FEATURES

CCL 1/2



Compression Class

Silicone beaded Top Band

Clima Fresh

Perfect Fit

Graduated Compression



Calf (Beige) M.R.P. Rs 3680/-pair
Thigh (TB-Beige) M.R.P. Rs 5688/-pair, Panty (Beige) M.R.P. Rs 6367/-pc

QUALITY CERTIFICATES

CE mark

Quality label (RAL)

Textile trust Quality seal

HOSY-Hohenstein Testing Device

Certificate
DIN EN ISO 9001:2008
DIN EN 13485

AVAILABLE MODELS

Calf

Thigh


Panty

Close toe

Open toe

MEASUREMENT

Sizes	cB Ankle measurement	cC Calf measurement	cG Thigh measurement	measurements for silicon topband
S	19 - 21 cm	28 - 34 cm	42 - 57 cm	47 - 54 cm
M	22 - 24 cm	32 - 38 cm	48 - 64 cm	52 - 60 cm
L	25 - 27 cm	36 - 42 cm	54 - 71cm	56 - 65 cm
XL	28 - 30 cm	40 - 46 cm	60 - 78 cm	62 - 72 cm
XXL	30 - 32 cm	42 - 50 cm	65 - 85 cm	65 - 76 cm



Correct graduation of Pressure depends on correct selection of the ankle size

(I) Emter. M, Compression Stocking management of patients suffering from Venous Disorder (II) Phlebology Net(5 Abstracts) (Stockings)



PUSHPANJALI medi India Pvt Ltd

MEDICAL COMPRESSION | ORTHOPEDICS | PHYSIOTHERAPY | REHABILITATION

