

mediven[®] ulcer kit/refill

medi

Effective compression therapy for the treatment of leg ulcer

Components :
2 mediven[®] ulcer
& 1 mediven[®] ulcer plus



Liner stocking mediven[®]ulcer

Therapeutically effective compression with 20 mmHg during the night.

Flat topband will stay in place without constriction

Soft material with antibacterial effect for pleasant wearing comfort and easy donning

Sewn-in instep area allows for optimum fit and reduces wrinkling when foot is bent

Special heel construction offers optimum pressure relief and colour marking for easy positioning

NEW

Outer stocking mediven[®]ulcer plus

Additional compression of 20 mmHg during the day.

Wide topband ensures optimum wearing comfort

Medium material strength for comfortable yet effective working pressure

Comfortable plus-heel and reinforced foot sole for good durability



2 x 20 mmHg
COMPRESSION CLASS

QUALITY CERTIFICATES

CE mark Quality label (RAL) Textile trust Quality seal HOSY-Hohenstein Testing Device DIN EN ISO 9001:2008 DIN EN 13485

QUALITY FEATURES

medi since- 1951 Soft topband Heel / instep region Antibacterial Easy donning Softness

HIGHLIGHTS

Double layer compression stocking with 40 (20+20) mmHg for the treatment of venous leg ulcers with material composition:

- mediven[®]ulcer: polyamide 73%, elasthane 27% (anti-microbial yarn)
- mediven[®]ulcer plus: polyamide 75%, elasthane 25%

(1) Pangratis.N, 1998, Chronic Venous ulcers- Conservativ Treatment-New trends

MEASUREMENT

Size	Leg measurements (in cm)			Colour mark
	Ankle cB	Calf cC	Below knee cD	
I	18-20	28-34	27-33	pink
II	20-22	30-37	29-36	blue
III	22-24	33-40	32-39	yellow
IV	24-26	35-43	34-42	red
V	26-28	37-46	36-45	grey
VI	28-30	39-49	38-48	green
VII	30-32	41-51	40-50	brown



PUSHPANJALI medi India Pvt Ltd
MEDICAL COMPRESSION | ORTHOPEDICS | PHYSIOTHERAPY | REHABILITATION



*All prices are subject to change without notice, and are not guaranteed. Please check during purchase / sale.